

# SEBRA<sup>®</sup>

## Hand held RF tube sealers

INTEGRATED DEVICES



# Model 2380 sealer



The battery operated Model 2380 hand held sealer is ideal for use in apheresis applications, for leukoreduction processes and anywhere mobility and electrical safety cannot be compromised. The 2380 hand held sealer comes with an 1105 sealing head, one NiCad battery, and a canvas carrying case.

- Certified for donor connected use
- Designed for mobile use
- 1500 seals per fully charged battery pack
- Two hour quick recharge
- Patented SMART™ Electronics automatically adjust for different size tubing; including whole blood and apheresis
- Tear seal feature makes segments that can be easily separated by hand
- Can be used by or near persons with pacemakers\*
- 20+ years proven field experience

# Model 2490 Sealite™ sealer

The battery operated Model 2490 Sealite tube sealer seals tubing used in apheresis collections and leukoreduction processing. The Sealite tube sealer is ideal for and anywhere mobility and electrical safety are important. The Sealite tube sealer includes power source with integrated battery, charger, carry case, coiled RF cable and Model 1105 Hand Held Sealing Head.

- Certified for donor connected use
- Tear seal feature makes segments that can be easily separated by hand
- 250 seals on mobile apheresis tubing per charge
- 325 seals on standard blood bag tubing per charge
- Can be used by or near persons with pacemakers\*



\*Results of testing on SEBRA Tube Sealers demonstrated that all SEBRA Tube Sealers tested, used in intended methods of application, comply with the MPE (maximum permissible exposure) limits specified for uncontrolled environments for energy density, electric field strength, and magnetic field strength under IEEE C95.1-1992, *Standard for Safety levels with Respect to Radio Frequency Electromagnetic Fields, 3kHz to 300GHz*

## Battery sealer comparison chart

Model Number	RoHS Compliant	Battery type	Type of tubing sealed	Seals per charge (on blood bag tubing)	Size (L x W x H)	Weight
2490	No	NiCad	blood bag, apheresis	325	Power Source 7.5 x 4.25 x 1.75 in. (20 x 11 x 5 cm) Sealing Head 8.5 x 1 in. (22 x 2.5 cm)	Power Source 2.4 lbs (1.1 kg) Sealing Head 0.6 lbs (.3 kg)
2380	No	NiCad	blood bag, apheresis, cell wash	1500	Power Source 12.5 x 4.8 x 1.7 in. (31.8 x 12.2 x 4.3 cm) Sealing Head 8.5 x 1 in. (22 x 2.5 cm)	Power Source 6 lbs (2.7 kg) Sealing Head 0.6 lbs (.3 kg)



Copyright © 2010-2011 Haemonetics Corporation.  
Haemonetics, Sealite, SMART and SEBRA are trademarks  
or registered trademarks of Haemonetics Corporation  
in the United States, other countries or both. All rights  
reserved. 05.2011 Germany. COL-COPY-000258(AB)