


# EIDorado Donor<sup>®</sup>

Comprehensive blood  
management system

Designed by blood bankers.  
For blood bankers.  
With your needs in mind.

SOFTWARE SOLUTIONS



The screenshot displays a donor profile form with the following fields and values:

- First: ADAM
- Middle: [Empty]
- Nick: [Empty]
- Preferred: [Empty]
- Birth Date: 8/19/1958
- Gender: M (Male)
- Marital: MARRIED (Married)
- Ethnicity: A (Asian)
- E-Mail: [Empty]
- Occupation: COMPUTER (Computer Analyst)
- Employer: [Empty]
- Language: U (Unknown)



# Comprehensive blood management system

The screenshot displays a software interface for blood management. At the top, it shows 'Component Manufacture Information' with fields for 'FE99902', '04210', and '10/15/2007 1:49 PM'. Below this, there are sections for 'Registration No.' (777777), 'Weight', 'Volume: 227.00', and 'Equipment ID'. A 'Draw Information' section includes 'Draw ID: 1410', 'Draw Location: BB2', 'Draw Date: 10/15/2007', 'Start Time: 1:49 PM', 'Stop Time: 1:50 PM', and 'Draw Status: C'. The interface also shows 'Blood Bank 2 - Alpha' and 'Consumed' status.

## Integrated modular solutions for

- System management
- Donor management
- Laboratory management
- Manufacturing
- Distribution

## System highlights

- Integrated modules<sup>1</sup> provide tailored solutions
- System scalability supports site and system growth
- Automated safety checks help control product release
- At-A-Glance-Bars allow easy access to critical data
- Multi-tabbed forms access all data with intuitive navigation at your fingertips
- Flexible component manufacturing processes increase efficiency in your lab
- Facility-specific SOP links can be tailored to your system
- ISBT 128 and ABC codabar automated blood labeling

## Technical platform

- Windows® operating system or Oracle® RDBMS
- Centralized security
- Auto-update of upgrades
- Crystal Reports®
- Centralized administration

The EIDorado suite of blood management products was developed using advanced tools and technologies to support the integrated data processing requirements of small centers to large complex organizations.

# Be Safe. Be Sure. Be Secure.

## Module highlights

### System management

- **Auditor**  
Track changes to data with secure time-stamped audit trails to independently record the date and time of operator entries and actions that create, modify or delete electronic records.
- **Advisor**  
Import and link customized SOP documents to forms, tabs, reports and messages for easy access to SOP content.
- **Security**  
Determine security access with multi-level sitedefined hierarchy allowing management to control access to critical data and set security paradigms for users and groups.
- **Table administration**  
Enjoy simple access to all table administration forms and user configuration of site business rules.
- **Central upgrade management**  
Automatically upgrade the application on all workstations from one central location.

### Donor management

- **Configurable pre-donation information**  
Gather information and record to consistently store data and help achieve accurate results.
- **Automatic blood loss and eligibility calculations**  
Apply configurable donor eligibility settings including automated blood loss tracking and view the status onscreen immediately to mitigate unnecessary draws.
- **Configurable donor identifier types**  
Store donor photos, identify donors using multiple methods and use a duplicate search mechanism to help preserve the integrity of donor identification.
- **Special collections**  
Automatically link donations from special collection requests to the patient and maintain that link throughout manufacturing, labeling and distribution.

### Manufacturing

- **Flexible component manufacturing processes**  
Increase efficiency in your lab, meet site specific operational needs, and enable staff to capture detailed manufacturing information.

### Laboratory management

- **Configurable testing**  
Streamline internal processes and save technician time using configurable testing which allows the system to automate test order generation based on intended use.
- **Site-defined reflex testing**  
Help prevent the manufacture of units until additional tests are completed with site-defined reflex testing.
- **Donor work queue**  
Save technician time with one-click access to appropriate testing functions and displays of pending test orders and tests.
- **Automatic safety checks**  
Performed throughout the system as part of the normal processing of donors, blood components and other blood products, to help prevent the release of unsuitable blood components.

### Inventory / distribution

- **Complete activity record**  
Track products and supplies more efficiently, as items are shipped or used in the system, and view the information easily in multi-tabbed profile forms.
- **Central intake**  
Receive an up-to-date inventory of available products at all sites, track products from the distribution location to customers or sites.
- **Configurable product options**  
Uniquely identify, inventory products, and create ISBT 128 or ABC codabar labeling at the product level assuring advanced traceability and organized distribution of products.



Our team embraces a company-wide Service 360® commitment to provide you with well-developed products, comprehensive education and professional resources, to enable you to utilize your products to their fullest. Our customer support analysts are always available to serve you.

## Key features and benefits:

### Safety features

- Identify duplicate donors
- Calculate donor eligibility
- Automate blood loss calculations
- Link comprehensive labeling and release rules to donor information, manufacturing and testing
- Track special collection requests

### Gain flexibility while meeting regulations

- Manufacture with flexibility
- Interchange ABC codabar and ISBT 128
- Access SOPs and help messages
- Perform routine safety checks automatically
- Tailor a system with integrated modules

### Increase efficiency

- Streamline the pre-donation process
- Access information quickly collect, track and store comprehensive donor and component information
- Use automated labeling and verification processes
- Manage inventory with multi-tabbed profile forms
- Access donor unit profiles with one-click feature

### United States

Haemonetics Corporation  
Software Solutions  
+1.800.996.3428  
[www.wyndgate.com](http://www.wyndgate.com)

### Europe | ME | Asia

INLOG S.A.S.  
A Haemonetics Company  
+33.478.66.53.53  
[www.inlog.com](http://www.inlog.com)



**Rx Only**

© 2010-2011 Haemonetics Corporation. Haemonetics, Eldorado, Eldorado Donor, Service 360, and Service 360 design are trademarks or registered trademarks of Haemonetics Corporation in the USA, other countries, or both. Crystal Reports is a registered trademark of Business Objects S.A. Oracle is a trademark of Oracle International Corporation. Windows is a trademark of Microsoft Corporation. 12.2011 Germany. COL-PP-000105-US(AB)